## Introducing



# Mona Lisa® Copper Intrauterine Devices

For the first time in Australia and New Zealand, Mona Lisa® NV has partnered with Basin Medical Pty Ltd to provide a new range of Copper IUD's. Germanmade Mona Lisa® Copper IUD's have been in production for over 30 years.

Mona Lisa Intrauterine Device	Width (mm)	Length (mm)	Maximum Insertion Time in the Body	Recommended Uterus Length
Mona Lisa Cu375	19.5	35.2	5 years	6 – 9 cm
Mona Lisa Cu375 SL	19.5	29.4	5 years	5 – 8 cm
Mona Lisa CuT 380A QL	31.9	35.9	10 years	6 – 9 cm

Barium Sulphate for radiological proofing.

Highest quality German made devices at lower prices.

MR Conditional.



#### Mona Lisa Cu375

• Copper surface: 375 mm<sup>2</sup>

In-situ-time: 5 years

One-handed insertion technique

#### Mona Lisa Cu375 SL

- Identical characteristics to the Mona Lisa Cu375
- Shorter stem length (29.4 mm)
- Ideally suited for a shorter uterus length.

#### Mona Lisa CuT 380A QL

Copper surface: 380 mm<sup>2</sup>

In-situ-time: 10 years

- IUD with the lowest Pearl Index
- Easy and quick loading of the IUD due to the new technology



Please contact Basin for more information.

# Mona Lisa® 🛕 Copper Intrauterine Devices

For the first time in Australia and New Zealand, Mona Lisa® NV has partnered with Basin Medical Pty Ltd to provide a new range of Copper IUD's. Germanmade Mona Lisa® Copper IUD's have been in production for over 30 years.

## Quality

Mona Lisa Copper IUD's are manufactured according to ASTM
Standards for MRI use, ISO 9001:2008, ISO 13485:2003, CAN/CSA ISO 13485:2003 and EN ISO 13485:2012





# Mona Lisa®

# Copper Intrauterine Devices

Mona Lisa® Copper IUD's, your MR Conditional alternative.

### What makes the Basin medical Mona Lisa Copper IUD's unique?

Basin Medical Copper IUD's are MR Conditional.

Non-clinical testing has demonstrated the Mona Lisa intrauterine device, represented by CuT 380 A QL as worst case, is MR conditional. A patient with this device can be safely scanned in an MR system meeting the certain conditions including a static magnetic field of 1.5 Tesla and 3 Tesla.



## Quality

Mona Lisa Copper IUD's are manufactured according to ASTM
Standards for MRI use, ISO 9001:2008, ISO 13485:2003, CAN/CSA ISO 13485:2003 and EN ISO 13485:2012

### Value

We pride ourselves on providing the highest quality German made Copper IUD's at lower prices. Contact us for a quotation.

Phone: 0429 130 945

## **Patient Safety**

More than 30 years of manufacturing expertise has seen Mona Lisa® Copper IUD's become one of the leading generic copper IUDs, worldwide.

For further information or to order please visit www.basinmedical.com.au

