


INTRODUCING...

Mona Lisa[®]

CuT 380A QL



FROM **BASIN**  **MEDICAL**


CuT 380A QL

This model is an enhanced version of the no longer produced “GyneT” with a new technology for quick and safe loading of the IUD.

The classic T-shaped IUD is the most common used IUD-model in the world, recommended even by the WHO.

With a copper surface of 380 mm² it shows the lowest pearl index among all copper IUDs.

The Mona Lisa CuT 380A QL ensures a secure fit in the uterus which produces a low expulsion rate due to optimal memory test results.



Unique features:

- ✓ The Mona Lisa CuT 380A has a copper surface of 380 mm² and guarantees safe contraception up to 10 years (recommendation of the World Health Organisation)
- ✓ The copper sleeves are completely embedded in the horizontal arms. This results in a very smooth and plane surface
- ✓ The longer insertion tube of the Mona Lisa CuT380A QL improves positioning of the IUD in patients with a longer uterine cavity
- ✓ Effortless and quick loading of the IUD because of the new insertion body
- ✓ The only Australian MR Conditional Copper IUD

BASIN  **MEDICAL**

For further information on these products visit the Basin Medical website www.basinmedical.com.au or email sales@basinmedical.com.au